



## Safty syringe with retractable needle

SIZE	0.5ML 1.0 ML 3ML/5ML/10ML
Features	1)automatically destroyed after single use
	2)after injection, its piston is automatically destroyed while the needle retreats plunger.
	3)once the plunger has been broken off, hence offering medical protection for both medical workers and patients,
Pakage	blister packing. Different package according to your requests.
	1)On the basis of traditional syringe, it has fuctions of auto-destroy structure and easy operation
	2)Can be destroyed automatically after single use

## RIBBEL® GERMAN MEDICAL DEVICES UG PLESTIA HEALTH CARE SERVICES

Tel: 03-7220577 Fax: 037219984